




BCAS2 Antibody

Product Code	CSB-PA446952
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75934
Immunogen	Synthesized peptide derived from C-terminal of Human BCAS2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BCAS2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Pre-mRNA-splicing factor SPF27; Spliceosome-associated protein SPF 27; Breast carcinoma-amplified sequence 2; DNA amplified in mammary carcinoma 1 protein; BCAS2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BCAS2

Image



Western blot analysis of extracts from HeLa cells, HepG2 cells, Huvec cells, A549 cells, K562 cells and Jurkat cells, using BCAS2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.