





BCAS4 Antibody

Product Code	CSB-PA225925
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDM0
Immunogen	Synthesized peptide derived from internal of Human BCAS4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BCAS4 protein.
Tested Applications	ELISA, WB, IHC; WB:1:500-1:3000, IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Breast carcinoma-amplified sequence 4; BCAS4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BCAS4

Image



Western blot analysis of extracts from K562 cells, using BCAS4 antibody.



Immunohistochemistry analysis of paraffin-embedded human tonsil tissue using BCAS4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.