




BAZ2A Antibody

Product Code	CSB-PA437543
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIF9
Immunogen	Synthesized peptide derived from internal of Human BAZ2A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total BAZ2A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	BA2A; BAZ2A; bromodomain adjacent to zinc finger domain 2A; hWALp3; KIAA0314
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAZ2A

Image



Western blot analysis of extracts from Jurkat cells, using BAZ2A antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.