




RCAN1 Antibody

Product Code	CSB-PA978873
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53805
Immunogen	Synthesized peptide derived from C-terminal of Human RCAN1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total RCAN1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADAPT78; calcipressin 1; CCP1; Down syndrome critical region protein 1; MCIP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RCAN1

Image



Western blot analysis of extracts from HeLa cells, using RCAN1 antibody.