




RCAN1 Antibody

Product Code	CSB-PA978873
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53805
Immunogen	Synthesized peptide derived from C-terminal of Human RCAN1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total RCAN1 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ADAPT78; calcipressin 1; CCP1; Down syndrome critical region protein 1; MCIP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RCAN1

Image



Western blot analysis of extracts from HeLa cells, using RCAN1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.