





CABP7 Antibody

Product Code	CSB-PA180621
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86V35
Immunogen	Synthesized peptide derived from internal of Human CABP7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CABP7 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Calcium-binding protein 7; CaBP7;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CABP7

Image



Western blot analysis of extracts from Jurkat cells, using CABP7 antibody.



Immunohistochemistry analysis of paraffin-embedded human heart tissue using CABP7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.