





CNN2 Antibody

Product Code	CSB-PA901238
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99439
Immunogen	Synthesized peptide derived from internal of Human CNN2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total CNN2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Calponin H2; smooth muscle; CLP2; CNN2; Neutral calponin
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNN2

Image



Western blot analysis of extracts from HUVEC cells, using CNN2 antibody.



Immunofluorescence analysis of A549 cells, using CNN2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.