




CALB2 Antibody

Product Code	CSB-PA077274
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22676
Immunogen	Synthesized peptide derived from internal of Human CALB2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CALB2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Calretinin; 29 kDa calbindin; CAB29; calretinin; CR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CALB2

Image



Western blot analysis of extracts from Jurkat cells, using CALB2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.