




ABCC3 Antibody

Product Code	CSB-PA105454
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15438
Immunogen	Synthesized peptide derived from internal of Human ABCC3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total ABCC3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Canalicular multispecific organic anion transporter 2; ATP-binding cassette sub-family C member 3; Multi-specific organic anion transporter D; MOAT-D; Multidrug resistance-associated protein 3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABCC3

Image



Western blot analysis of extracts from HeLa cells, A549 cells, Jurkat cells, HuvEc cells, K562 cells and HepG2 cells, using ABCC3 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.