




CLIP1 Antibody

Product Code	CSB-PA183794
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30622
Immunogen	Synthesized peptide derived from internal of Human CLIP1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total CLIP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CLIP-170; CLIP1; Cytoplasmic linker protein-170 alpha-2; Reed- Sternberg intermediate filament associated protein; REST
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLIP1

Image



Western blot analysis of extracts from MCF-7 cells, RAW264.7 cells, HepG2 cells, COLO cells and HUVEC cells, using CLIP1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.