




CLIP4 Antibody

Product Code	CSB-PA043960
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N3C7
Immunogen	Synthesized peptide derived from internal of Human CLIP4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CLIP4 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Canalicular multispecific organic anion transporter 1; ATP-binding cassette sub-family C member 2; Canalicular multidrug resistance protein; Multidrug resistance-associated protein 2; ABCC2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLIP4

Image



Western blot analysis of extracts from HepG2 cells and Jurkat cells, using CLIP4 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.