





CHST10 Antibody

Product Code	CSB-PA877976
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43529
Immunogen	Synthesized peptide derived from internal of Human CHST10.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CHST10 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Carbohydrate sulfotransferase 10; EC 2.8.2.-; HNK-1 sulfotransferase; HNK1ST; HNK-1ST
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CHST10

Image



Western blot analysis of extracts from HuvEc cells, using CHST10 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using CHST10 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.