




HTR3A Antibody

Product Code	CSB-PA786914
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46098
Immunogen	Synthesized peptide derived from internal of Human 5-HT-3A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total 5-HT-3A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	5-HT-3; 5-HT3R; 5-hydroxytryptamine 3 receptor; 5HT3R; HTR3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HTR3A

Image



Western blot analysis of extracts from HT-29 cells, using 5-HT-3A antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.