





CA1 Antibody

Product Code	CSB-PA545638
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00915
Immunogen	Synthesized peptide derived from internal of Human CA1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CA1 protein.
Tested Applications	ELISA, WB, IHC; WB:1:500-1:3000, IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Carbonic anhydrase 1; EC 4.2.1.1; Carbonic anhydrase I; Carbonate dehydratase I; CA-I
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CA1

Image



Western blot analysis of extracts from K562 cells, using CA1 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using CA1 antibody.