




CA14 Antibody

Product Code	CSB-PA962891
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULX7
Immunogen	Synthesized peptide derived from internal of Human CA14.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CA14 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Carbonic anhydrase 14; EC 4.2.1.1; Carbonic anhydrase XIV; Carbonate dehydratase XIV; CA-XIV
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CA14

Image



Western blot analysis of extracts from 293 cells, COLO cells, and Jurkat cells, using CA14 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.