





# CA5B Antibody

<b>Product Code</b>	CSB-PA039070
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y2D0
<b>Immunogen</b>	Synthesized peptide derived from internal of Human CA5B.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous levels of total CA5B protein.
<b>Tested Applications</b>	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:200,IF:1:100-1:200
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Carbonic anhydrase 5B; mitochondrial [Precursor]; EC 4.2.1.1; Carbonic anhydrase VB; CA-VB
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CA5B

## Image



Western blot analysis of extracts from NIH/3T3 cells, using CA5B antibody.



Immunofluorescence analysis of MCF-7 cells, using CA5B antibody.

### Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.