





CA5B Antibody

Product Code	CSB-PA039070
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2D0
Immunogen	Synthesized peptide derived from internal of Human CA5B.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CA5B protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:200,IF:1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Carbonic anhydrase 5B; mitochondrial [Precursor]; EC 4.2.1.1; Carbonic anhydrase VB; CA-VB
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CA5B

Image



Western blot analysis of extracts from NIH/3T3 cells, using CA5B antibody.



Immunofluorescence analysis of MCF-7 cells,
using CA5B antibody.