




CA6 Antibody

Product Code	CSB-PA167607
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23280
Immunogen	Synthesized peptide derived from internal of Human CA6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CA6 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Carbonic anhydrase 6; EC=4.2.1.1; Carbonic anhydrase VI; CA-VI; Carbonate dehydratase VI
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CA6

Image



Western blot analysis of extracts from LOVO cells and HepG2 cells, using CA6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.