





CBR3 Antibody

Product Code	CSB-PA669024
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75828
Immunogen	Synthesized peptide derived from internal of Human CBR3.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CBR3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	carbonyl reductase 3; EC 1.1.1.184; NADPH-dependent carbonyl reductase 3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CBR3

Image



Western blot analysis of extracts from MCF7 cells and Jurkat cells, using CBR3 antibody.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue, using CBR3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.