




CPA5 Antibody

Product Code	CSB-PA011345
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Swiss-Prot: Q8WXQ8 NCBI Gene ID: 93979
Immunogen	Synthesized peptide derived from internal of Human CPA5.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CPA5 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Carboxypeptidase A5;
Immunogen Species	Homo sapiens (Human)
Target Names	CPA5

Image



Western blot analysis of extracts from 293 cells, using CPA5 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.