




CEACAM21 Antibody

Product Code	CSB-PA942626
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3KPI0
Immunogen	Synthesized peptide derived from internal of Human CEACAM21.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CEACAM21 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Carcinoembryonic antigen-related cell adhesion molecule 21;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEACAM21

Image



Western blot analysis of extracts from K562 cells, using CEACAM21 antibody.