



CEP152 Antibody

Product Code	CSB-PA042364
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94986
Immunogen	Synthesized peptide derived from internal of Human CEP152.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total CEP152 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	centrosomal protein 152 kDa; centrosomal protein 152kDa; centrosomal protein of 152 kDa;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEP152

Image



Immunohistochemistry analysis of paraffin-embedded human lymph node tissue using CEP152 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.