




CEP76 Antibody

Product Code	CSB-PA904419
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAP6
Immunogen	Synthesized peptide derived from internal of Human CEP76.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CEP76 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C18orf9; centrosomal protein 76kDa; HsT1705;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEP76

Image



Western blot analysis of extracts from 3T3 cells, using CEP76 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.