





CLN6 Antibody

Product Code	CSB-PA167979
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWW5
Immunogen	Synthesized peptide derived from internal of Human CLN6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CLN6 protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ceroid-lipofuscinosis neuronal protein 6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLN6

Image



Western blot analysis of extracts from HeLa cells, using CLN6 antibody.



Immunohistochemistry analysis of paraffin-embedded human cervix tissue using CLN6 antibody.