





CLCN4 Antibody

Product Code	CSB-PA968321
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51793
Immunogen	Synthesized peptide derived from internal of Human CLCN4.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total CLCN4 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:200,IF:1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	chloride channel 4; CIC-4; CIC-4A; CLC4; MGC163150
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLCN4

Image



Western blot analysis of extracts from MCF-7 cells, using CLCN4 antibody.



Immunofluorescence analysis of HUVEC cells, using CLCN4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.