




CLCN7 Antibody

Product Code	CSB-PA212493
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51798
Immunogen	Synthesized peptide derived from N-terminal of Human CLCN7.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CLCN7 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Chloride channel protein 7; CIC-7;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CLCN7

Image



Western blot analysis of extracts from A549 cells and LOVO cells, using CLCN7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.