




CDYL2 Antibody

Product Code	CSB-PA019829
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N8U2
Immunogen	Synthesized peptide derived from N-terminal of HumanCDYL2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CDYL2 protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Chromodomain Y-like protein 2; CDY-like 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDYL2

Image



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using CDYL2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.