




CCDC102A Antibody

Product Code	CSB-PA065026
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A19
Immunogen	Synthesized peptide derived from internal of Human CCDC102A.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CCDC102A protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C102A; coiled-coil domain-containing 102A;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCDC102A

Image



Western blot analysis of extracts from LOVO cells, using CCDC102A antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.