




CCDC102B Antibody

Product Code	CSB-PA190186
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68D86
Immunogen	Synthesized peptide derived from internal of Human CCDC102B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CCDC102B protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 102B;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCDC102B

Image



Western blot analysis of extracts from HuvEc cells and COLO cells, using CCDC102B antibody.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using CCDC102B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.