




C1QB Antibody

Product Code	CSB-PA197114
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02746
Immunogen	Synthesized peptide derived from internal of Human C1QB.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total C1QB protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:200,IF:1:100-1:200
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	complement C1q subcomponent B; complement component 1; q subcomponent; B;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	C1QB

Image



Western blot analysis of extracts from Jurkat cells, using C1QB antibody.



Immunofluorescence analysis of A549 cells, using C1QB antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.