




CYP2A6 Antibody

Product Code	CSB-PA909765
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11509
Immunogen	Synthesized peptide derived from N-terminal of Human Cytochrome P450 2A6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 2A6 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450 2A6; EC 1.14.14.1; CYP11A6; Coumarin 7-hydroxylase; P450 IIA3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2A6

Image



Western blot analysis of extracts from Jurkat cells, using Cytochrome P450 2A6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.