




CYP2A6 Antibody

Product Code	CSB-PA909765
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11509
Immunogen	Synthesized peptide derived from N-terminal of Human Cytochrome P450 2A6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total Cytochrome P450 2A6 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cytochrome P450 2A6; EC 1.14.14.1; CYP11A6; Coumarin 7-hydroxylase; P450 IIA3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CYP2A6

Image



Western blot analysis of extracts from Jurkat cells, using Cytochrome P450 2A6 antibody.