




CNTN4 Antibody

Product Code	CSB-PA068627
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWW2
Immunogen	Synthesized peptide derived from internal of Human CNTN4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total CNTN4 protein.
Tested Applications	ELISA, IF; IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Contactin-4; Brain-derived immunoglobulin superfamily protein 2; BIG-2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNTN4

Image



Immunofluorescence analysis of A549 cells, using CNTN4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.