




CNTN6 Antibody

Product Code	CSB-PA973774
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQ52
Immunogen	Synthesized peptide derived from internal of Human CNTN6.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CNTN6 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	contactin 6; contactin-6; NB3; neural adhesion molecule; neural recognition molecule NB-3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNTN6

Image



Western blot analysis of extracts from Jurkat cells and A549 cells, using CNTN6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.