





HSD11B1 Antibody

Product Code	CSB-PA213517
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28845
Immunogen	Synthesized peptide derived from N-terminal of Human HSD11B1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total HSD11B1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Corticosteroid 11-beta-dehydrogenase isozyme 1; 11-DH; 11-beta-hydroxysteroid;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSD11B1

Image



Western blot analysis of extracts from HuvEc cells, using HSD11B1 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using HSD11B1 antibody.