




CXADR Antibody

Product Code	CSB-PA830526
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78310
Immunogen	Synthesized peptide derived from internal of Human CXADR.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CXADR protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CAR; Coxsackievirus and adenovirus receptor precursor; Coxsackievirus B-adenovirus receptor; CVB3 binding protein; CXAR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CXADR

Image



Western blot analysis of extracts from COLO cells and HepG2 cells, using CXADR antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.