



CNGA2 Antibody

Product Code	CSB-PA139476
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16280
Immunogen	Synthesized peptide derived from internal of Human CNGA2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total CNGA2 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CNCA; CNCA1; CNCG2; CNG channel 2; CNG-2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CNGA2

Image



Immunofluorescence analysis of A549 cells, using CNGA2 antibody.