



CSTL1 Antibody

Product Code	CSB-PA030996
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H114
Immunogen	Synthesized peptide derived from internal of Human CSTL1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CSTL1 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cystatin-like 1; RCET11;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CSTL1

Image



Immunofluorescence analysis of A549 cells, using CSTL1 antibody.