




CST1 Antibody

Product Code	CSB-PA963535
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01037
Immunogen	Synthesized peptide derived from internal of Human CST1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CST1 protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	cystain-SA-I; cystatin-1; cystatin-SN; CYTN; salivary cystatin-SA-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CST1

Image



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using CST1 antibody.



Immunofluorescence analysis of MCF-7 cells,
using CST1 antibody.