




NT5C1B Antibody

Product Code	CSB-PA969848
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P26
Immunogen	Synthesized peptide derived from internal of Human NT5C1B.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total NT5C1B protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	5 nucleotidase; cytosolic IB; 5NT1B; AIRP; EC 3.1.3.5
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NT5C1B

Image



Immunofluorescence analysis of A549 cells, using NT5C1B antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.