




NT5C1B Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA969848 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96P26 |
| Immunogen | Synthesized peptide derived from internal of Human NT5C1B. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Specificity | The antibody detects endogenous levels of total NT5C1B protein. |
| Tested Applications | ELISA,IF;IF:1:100-1:500 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | 5 nucleotidase; cytosolic IB; 5NT1B; AIRP; EC 3.1.3.5 |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | NT5C1B |

Image



Immunofluorescence analysis of A549 cells, using NT5C1B antibody.