




APPL1 Antibody

Product Code	CSB-PA246914
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKG1
Immunogen	Synthesized peptide derived from internal of Human APPL1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total APPL1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Adaptor protein APPL; Adaptor protein containing pH domain; PTB domain and leucine zipper motif; DIP13 alpha; DP13A
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	APPL1

Image



Western blot analysis of extracts from A549 cells, using APPL1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.