




DHRS2 Antibody

Product Code	CSB-PA273664
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13268
Immunogen	Synthesized peptide derived from internal of Human DHRS2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total DHRS2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	dehydrogenase/reductase (SDR family) member 2; dehydrogenase/reductase SDR family member 2; HEP27; protein D; SDR25C1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DHRS2

Image



Western blot analysis of extracts from COLO cells, using DHRS2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.