




POLG2 Antibody

Product Code	CSB-PA826673
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHN1
Immunogen	Synthesized peptide derived from internal of Human POLG2.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total POLG2 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA polymerase subunit gamma-2; mitochondrial; EC 2.7.7.7; Mitochondrial DNA polymerase accessory subunit; PolG-beta
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLG2

Image



Immunohistochemistry analysis of paraffin-embedded human cervix carcinoma tissue using POLG2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.