




PRIM1 Antibody

Product Code	CSB-PA972371
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49642
Immunogen	Synthesized peptide derived from C-terminal of Human PRIM1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PRIM1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA primase 1; DNA primase 49 kDa subunit; DNA primase small subunit; EC 2.7.7.-; PRI1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PRIM1

Image



Western blot analysis of extracts from HepG2 cells, using PRIM1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.