




POLR1B Antibody

Product Code	CSB-PA799360
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H9Y6
Immunogen	Synthesized peptide derived from C-terminal of Human POLR1B.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total POLR1Bprotein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA-directed RNA polymerase I 135 kDa polypeptide; DNA-directed RNA polymerase I 135kDa polypeptide; DNA-directed RNA polymerase I subunit RPA2; FLJ10816; FLJ21921
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	POLR1B

Image



Western blot analysis of extracts from Jurkat/NIH-3T3 cells, using POLR1B antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.