



MARCH4 Antibody

Product Code	CSB-PA973628
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2E8
Immunogen	Synthesized peptide derived from internal of Human MARCH4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total MARCH4 protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase MARCH4 ; EC 6.3.2.-; Membrane-associated RING finger protein 4; Membrane-associated RING-CH protein IV; MARCH-IV
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MARCH4