




RAD18 Antibody

Product Code	CSB-PA274794
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS91
Immunogen	Synthesized peptide derived from internal of Human RAD18.
Raised In	Rabbit
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous levels of total RAD18 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RAD18; hHR18; hRAD18; postreplication repair protein hRAD18p; Postreplication repair protein RAD18
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RAD18

Image



Western blot analysis of extracts from LOVO cells, HepG2 cells and K562 cells, using RAD18 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.