




RNF138 Antibody

Product Code	CSB-PA994526
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVD3
Immunogen	Synthesized peptide derived from internal of Human RNF138.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total RNF138 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF138; EC 6.3.2.-; RING finger protein 138; Nemo-like kinase-associated RING finger protein; NLK-associated RING finger protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	RNF138

Image



Western blot analysis of extracts from MCF-7 cells, using RNF138 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.