





EFEMP2 Antibody

Product Code	CSB-PA103283
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95967
Immunogen	Synthesized peptide derived from internal of Human EFEMP2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total EFEMP2 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EGF-containing fibulin-like extracellular matrix protein 2; FBLN4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EFEMP2

Image



Western blot analysis of extracts from K562 cells and 3T3 cells, using EFEMP2 antibody.



Immunohistochemistry analysis of paraffin-embedded human placenta tissue using EFEMP2 antibody.



Immunofluorescence analysis of A549 cells, using EFEMP2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.